



*Trans, The Highest Quality
Products For Life Sciences*

L-Glutamine (100×)

Cat. No. FG201

Storage: at -20°C for one year

Description

L-Glutamine (100×) contains 200 mM L-Glutamine. This product has been filter-sterilized. It can be directly used for cell culture based on a dilution ratio at 1:100.

Note

To avoid repeated freezing and thawing, aliquot after first thawing; store at 2-8°C for no more than two weeks.

FOR RESEARCH USE ONLY

Website www.transgenbiotech.com

E-mail trade@transgen.com.cn

Customer Service 0086-10-57815027

Phone 0086-10-57815087