



Product Information Sheet

M561 Murashige & Skoog (MS) Modified Basal Salt Mixture (1/2x Nitrates)

Properties

Form:	Fine to Fluffy Powder
Appearance:	White to Yellow Powder
Application:	Plant Tissue Culture
Solubility:	Water
Typical Working Concentration:	2.56 g/L
Storage Temp:	2-6°
Storage Temp of Stock Solution:	Preparation of concentrated solutions is not recommended as insoluble precipitates may form.
Other Notes:	Contains the macro- and micronutrients as described by Murashige and Skoog (1962) with the following exceptions: ½ Ammonium Nitrate and ½ Potassium Nitrate

Formula (mg/L)

Ammonium Nitrate	825
Boric Acid	6.2
Calcium Chloride, Anhydrous	332.2
Cobalt Chloride•6H ₂ O	0.025
Cupric Sulfate•5H ₂ O	0.025
Na ₂ EDTA•2H ₂ O	37.26
Ferrous Sulfate•7H ₂ O	27.8
Magnesium Sulfate, Anhydrous	180.7

Manganese Sulfate•H ₂ O	16.9
Molybdic Acid (Sodium Salt)•2H ₂ O	0.25
Potassium Iodide	0.83
Potassium Nitrate	950
Potassium Phosphate, Monobasic	170
Zinc Sulfate•7H ₂ O	8.6

Application Notes

Plant Tissue Culture Tested

References

Murashige, T and F Skoog. 1962. A revised medium for rapid growth and bioassays with tobacco tissue cultures. *Physiol. Plant.* 15: 473-497.

Revised 3/2007

PhytoTechnology Laboratories, LLC

P.O. Box 12205; Shawnee Mission, KS 66282-2205
Phone: 1-888-749-8682 or 913-341-5343; Fax: 1-888-449-8682 or 913-341-5442
Web Site: www.phytotechlab.com © 2007 PhytoTechnology Laboratories, LLC