



*Trans, The Highest Quality
Products For Life Sciences*

6×Protein Loading Buffer

Cat. No. DL101

Storage: at -20°C for one year

Description

This product is used as Loading Buffer for preparing protein sample in SDS-PAGE. Prior to loading, add appropriate volume of 6×Protein Loading Buffer to protein sample to make it working concentration at 1×.

FOR RESEARCH USE ONLY

Website www.transgenbiotech.com

E-mail trade@transgen.com.cn

Customer Service 0086-10-57815027

Phone 0086-10-57815087