



*Trans, The Highest Quality  
Products For Life Sciences*

## **Trypsin (0.25%, with EDTA and Phenol Red)**

Cat. No. FG301

Storage: at -20°C for 18 months

### **Description**

Trypsin contains porcine trypsin (0.25%) , EDTA, and phenol red. It does not contain calcium and magnesium ion. Trysin/EDTA solution has been filter-sterilized and it can be used for cell dissociation.

### **Note**

Packing after thawing; Store at 2-8°C not more than two weeks.

**FOR RESEARCH USE ONLY**