



*Trans, The Highest Quality
Products For Life Sciences*

***ProteinFind™* Anti-GAPDH Mouse Monoclonal Antibody**

Cat. No. HC301

Concentration: 1 mg/ml

Storage: at 2-8°C for one month, at -20°C for one year

Description

GAPDH (glyceraldehyde 3-phosphate dehydrogenase) is an essential enzyme that catalyzes the glycolysis, and has a high level of expression in most tissues and cells. It is commonly used as expression control. *ProteinFind™* Anti-GAPDH Mouse Monoclonal Antibody is purified mouse monoclonal IgG1 antibody that detects recombinant proteins containing GAPDH peptide sequence from a variety of species including human, mouse and rabbit.

Suggested Dilution

- Western Blot: 1:500-10000 dilution.
- ELISA: 1:500-5000 dilution.
- IF: 1:100-500 dilution.

Notes

- Aliquot prior to use, and avoid repeated freezing and thawing.
- Diluted antibody should be used within 12 hours, and avoid storing diluted antibody.

FOR RESEARCH USE ONLY

WeThink Biotech INC. (www.wethink-biotech.com)

wethinkbio@gmail.com (order), service.wethinkbio@gmail.com (support)