



*Trans, The Highest Quality
Products For Life Sciences*

***ProteinFind™* Goat Anti-Rabbit IgG(H+L), HRP Conjugate**

Cat. No. HS101

Concentration: 1 mg/ml

Storage: at 2-8°C for one month; at -20°C for one year

Description

Affinity purified *ProteinFind™* Goat Anti-Rabbit IgG(H+L) Antibody is conjugated to horseradish peroxidase. The product is used as a secondary antibody in western blot application.

Suggested Dilution

- Western Blot: 1:500-10000 dilution.
- ELISA: 1:500-5000 dilution.

Notes

- Aliquot prior to use, and avoid repeated freezing and thawing.
- Diluted antibody should be used within 12 hours, and avoid storing diluted antibody.
- Sodium azide is an inhibitor of HRP, thus avoid using sodium azide as a preservative in antibody dilution.

FOR RESEARCH USE ONLY

WeThink Biotech INC. (www.wethink-biotech.com)

wethinkbio@gmail.com (order), service.wethinkbio@gmail.com (support)