



*Trans, The Highest Quality  
Products For Life Sciences*

## ***ProteinFind*<sup>™</sup> Goat Anti-Rabbit IgG(H+L), PE Conjugate**

**Cat. No.** HS121

**Concentration:** 0.4 mg/ml

**Storage:** at 2-8°C in dark for one year

### **Description**

PE (phycoerythrin) is a natural fluorescent dye extracted from the red algae. PE has a maximum absorption wavelength of 488 nm and the greatest emission wavelength of 575 nm, with the orange red fluorescence. Affinity purified *ProteinFind*<sup>™</sup> Goat Anti-Rabbit IgG(H+L) Antibody is conjugated to PE. The product is used as a secondary antibody in application of immunofluorescence and flow cytometry.

### **Suggested Dilution**

- IF: 1:100-500 dilution.
- FCM: 1:100-500 dilution.

### **Notes**

- Aliquot prior to use, and store in dark.
- Diluted antibody should be used within 12 hours, and avoid storing diluted antibody.

**FOR RESEARCH USE ONLY**

WeThink Biotech INC. ([www.wethink-biotech.com](http://www.wethink-biotech.com))

[wethinkbio@gmail.com](mailto:wethinkbio@gmail.com) (order), [service.wethinkbio@gmail.com](mailto:service.wethinkbio@gmail.com) (support)