



*Trans, The Highest Quality
Products For Life Sciences*

***ProteinFind™* Anti-c-Myc Mouse Monoclonal Antibody**

Cat. No. HT101

Concentration: 1 mg/ml

Storage: at 2-8°C for one month, at -20°C for one year

Description

ProteinFind™ Anti-c-Myc Mouse Monoclonal Antibody is purified monoclonal antibody that detects recombinant proteins containing the c-Myc (EQKLISEEDL) amino acid sequence.

Suggested Dilution

- Western Blot: 1:500-10000 dilution.
- ELISA: 1:500-5000 dilution.
- IF: 1:100-500 dilution.
- IP: 1:100-500 dilution.

Notes

- Aliquot prior to use, and avoid repeated freezing and thawing.
- Diluted antibody should be used within 12 hours, and avoid storing diluted antibody.

FOR RESEARCH USE ONLY

WeThink Biotech INC. (www.wethink-biotech.com)

wethinkbio@gmail.com (order), service.wethinkbio@gmail.com (support)