



*Trans, The Highest Quality  
Products For Life Sciences*

## ***ProteinFind™ Anti-V5 Mouse Monoclonal Antibody***

Cat. No. HT401

Concentration: 1 mg/ml

Storage: at 2-8°C for one month, at -20°C for one year

### **Description**

*ProteinFind™* Anti-V5 Mouse Monoclonal Antibody is purified monoclonal antibody that detects recombinant proteins containing the V5 (CGKPIP NPLLGLDST) amino acid sequence.

### **Suggested Dilution**

- Western Blot: 1:500-10000 dilution.
- ELISA: 1:500-5000 dilution.
- IF: 1:100-500 dilution.
- IP: 1:100-500 dilution.

### **Notes**

- Aliquot prior to use, and avoid repeated freezing and thawing.
- Diluted antibody should be used within 12 hours, and avoid storing diluted antibody.

**FOR RESEARCH USE ONLY**

WeThink Biotech INC. ([www.wethink-biotech.com](http://www.wethink-biotech.com))

[wethinkbio@gmail.com](mailto:wethinkbio@gmail.com) (order), [service.wethinkbio@gmail.com](mailto:service.wethinkbio@gmail.com) (support)