



*Trans, The Highest Quality
Products For Life Sciences*

Penicillin-Streptomycin (100×)

Cat. No. FG101

Storage: at -20°C for one year

Description

Penicillin-Streptomycin (100×) contains 10 kU/ml of penicillin and 10 mg/ml of streptomycin. This product has been filter-sterilized. It can be directly used for cell culture based on a dilution ratio at 1:100.

Note

To avoid repeated freezing and thawing, aliquot after first thawing; store at 2-8°C for no more than two weeks.

FOR RESEARCH USE ONLY