



Digital fluorescence microscope

For RFP/mCherry

For DsRed

For CFP

For GFP



Fluorescent module for microscope

GFP/DsRed signal could be detected with Fluorescent Module.

GFP/DsRed signal could be detected with Fluorescent Module.



ZGENEBIO Stereo Microscope

USD\$1650

- **Lens:** binocular eyepiece with 45 degree heeled, spin 360 degree, and parallax of lens could be regulated for 52~75mm. The focusing knob could regulate the height.
- **Zoom:** 0.65X-4.5X
- **Amplification factor:** 10x Eyepiece, 6.5x-45x20x, and 13X-90X
Highest for 225x (diopter lens additionally purchasing)
- Suitable for anatomy, oocyte microinjection for fish, and drosophila
- Able to be connected with the digital image system
- **Transillumination base**
- **LED white light emission**
- **International design**



Transillumination base



Homepage of ZGene Biotech Inc. products
<http://zgb-microinjector.weebly.com/>